



Access for the patient is easy, (s)he can sit down easily, without being hindered by the footrest which is completely tucked away under the seat. The patient can also help him(her)self by resting on the armrests. The stability of the appliance is ensured by using the brakes on the wheels.

Moving a patient to the chair is made easy thanks to the thoughtful design: low profile dimensions behind the backrest which allow an ideal position for nursing staff and armrests collapsible to 180°.

The use of the footrest is very easy, it pulls out or pushes back in a single movement.

Mobility: the "HOMER" weighing chair equipped with 4 wheels (+ brakes) is easy to move. It also offers the advantage of being able to move with a patient sitting on the scales. The reliability, accuracy and simplicity of use of the latest electronic equipment serves patients and nursing staff.

It only takes three seconds to obtain accurate weighing results, the end of the weighing is signalled by a "beep".

The appliance is equipped with the following standard functions:

- RESET: restarts the system, returning to zero even if a load is present (enables the net weight to be determined).

- HOLD: maintains the weight display, during the whole duration of the weighing, from the time that the display is stabilised.

- REFRESH: restarts the weighing (HOT) without having to lift the patient from the scales and/or perform a stabilising of the display (anti-trembling).

- MEM: displays the result of a weighing which is stored in the memory.

The power supply can be chosen from the mains electric network (standard lead) or via an integrated battery (10 or 14 hours of battery life) to carry out up to several thousand weighings before recharging. A high performance control system and limit to the battery charge allow the appliance to be left plugged into the mains without any special precautions. Switching off is controlled manually or is automatic when running on battery.

## Technical Features

<b>Capacity / Accuracy :</b>	200kg / 100g
<b>Structure of the main chassis :</b>	reinforced thermoset resin
<b>Standard Functions :</b>	RESET (tare), HOLD, REFRESH, MEM Reset to zero (automatic) Anti-trembling Restart the current weighing Save the weighing Automatic turn off
<b>Electronic indicator :</b>	Behind the backrest
<b>Display :</b>	Digital LCD, easily readable
<b>Power supply (integrated) :</b>	100-240 Vac / 50-60 Hz / 200 mA
<b>Battery and charger :</b>	Built in (battery life 40h)
<b>Mobility :</b>	4 castors + brakes
<b>Armrests :</b>	Individually collapsible to 180°
<b>Footrest :</b>	Collapsible (17cm)
<b>Other well developed points :</b>	Hygienic – Ergonomic – Safe
<b>Back seat :</b>	40 x 46 cm
<b>Weighing surface/height :</b>	44 x 43 cm / 55 cm
<b>Total dimensions :</b>	59 x 55 x 100 cm
<b>Weight :</b>	22kg
<b>T° / HR accuracy guarantee: :</b>	from +5 °C to +40°C / <85%HR
<b>Manufacturers guarantee :</b>	2 years

## Options

<b>Capacity / Accuracy :</b>	200kg / 50g
------------------------------	-------------